



exocad ChairsideDB 3.1Rijeka

Exocam DWX-43W v25

Quick Guide

Introduction to exocam with SUM3D for DWX-43W Workflow

Workflow for DWX-43W Users

Workflow

Import the scan data of the prepared abutment teeth. Design the crown using exocad. Then create the machining data using exocam. The output machining data is processed through VPanel and machined using the DWX-43W.

Step1

Import intraoral scan data

Step2

exocad: Design Crown


Step3

exocam: CAM process

Step4

Milling

Workflow



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

Usage Conditions

[Requirements] In order to use this integration, the following requirements must be fulfilled.

Q 1 : Is VPanel for DWX installed?

No

Install VPanel for DWX.



Yes

Q 2 : Is exocad ChairsideCAD 3.1 Rijeka or later installed?

No

Update/install exocad Chairside CAD 3.1 Rijeka or later.



Yes

Q3 : Is DGSHAPE CAM V25.1.0 or later installed?

No

Update/install DGSHAPE CAM V25.1.0 or later.



Yes



Usage requirements are met. Please follow the next pages to continue the setup.



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

Preliminary Setup

*For installation settings, please contact the dealer.

[Installation EXOCAM]




Exocam DWX-43W v25

Unzip the file [EXOSUM3D_DWX-43W-HQ-2025_05_25.zip].

Install EXOSUM3D_DWX-43W-HQ-2025_05_25¥Setup.exe by right-clicking the file and selecting "Run as administrator".

Run EXOSUM3D_DWX-43W-HQ-2025_05_25/Setup.exe as Administrator to install the software.

Once the installation is complete, a shortcut icon named [Exocam DWX-43W v25] will be created on the desktop.



*This version of Exocam uses the SUM3D engine and is a CAM software that integrates seamlessly with exocad's CAD system, compatible with the DWX-43W.



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

Preliminary Setup

*For installation settings, please contact the dealer.

[Installation EXOCAD]

exocad ChairsideDB 3.1Rijeka

Regarding the Exocad ChairsideCAD installer, please obtain the “exocad-ChairsideCAD for DWX43W” installer directly from exocad GmbH.

To launch Exocad,
for example:

navigate to CHAIRSIDE EXOCAD_Roland DWX-43W_2025¥exocad-ChairsideCAD3.1-2023-09-12_DWX43W¥exocad-ChairsideCAD3.1-2023-09-12_DWX43W¥ChairsideDB¥bin¥

and double-click **DentalDB.exe**.

It is recommended to create a shortcut for DentalDB.exe and place it on the desktop.



WORKFLOW

IMPORT

DESIGN


CAM

MILLING

Preliminary Setup

*For LAN settings, please contact the dealer.

[LAN Setup]



(1) Connect the PC and the DWX-43W using a USB cable.



(2) Launch VPanel and select the DWX-43W you want to use.




(3) Click on [Machine Settings].



(4) In the LAN settings section, click [Add Connection] .



(5) Click [Next] .



(6) Select "Set up and connect the LAN".



(7) Confirm the IP address and click [Next] .



⑧ Click [Complete].



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

Step1 Import the intraoral scan data



WORKFLOW


IMPORT

DESIGN

CAM

MILLING

Step1 Import the intraoral scan data




Project

New Load Save Duplicate

Type in the patient name or select an existing patient. Hide names

2025/04/25 13:31 Due

Patient
DGS SHAPE



← Tooth 45 Material configuration (local): Default

Crowns and bridges

Anatomic crown Anatomic pontic

Inlays, onlays and veneers

Inlay/Onlay Veneer

Residual dentition

Antagonist Adjacent tooth Omit in bridge

More...

① First, register the patient.
Enter the patient name.

② Click on the tooth.

③ Select **Anatomic crown** (Anatomic crown).



WORKFLOW



IMPORT

DESIGN

CAM

MILLING

Step1 Import the intraoral scan data



- ④ Select the Material.
The materials selected here will be reflected in the materials of exocam, but you can change the materials in exocam.
 - ⑤ Set the 'Gap width of cement' and click the  (OK)



WORKFLOW


IMPORT

DESIGN

CAM

MILLING

Step1 Import the intraoral scan data



⑥ Return to the initial registration screen and click the  (Save) button at the top left of the screen. Next, click the  (Design) button at the top right.



WORKFLOW

IMPORT

DESIGN


CAM

MILLING

Step1 Import the intraoral scan data



⑦ Import the lower jaw scan data.



⑧ Import the upper jaw scan data.



WORKFLOW

IMPORT


DESIGN

CAM

MILLING

Step2 Design crown




Step2 Design crown



⑨ Adjust the orientation of the scan data.

Orientate the model data so that the occlusal view is seen.

If correct, click the **NEXT →** (Next) button to proceed.



⑩ In **Detect** (Detect) function, click the margin line to start automatic detection. If automatic detection fails, delete the margin line with **CLEAR** (CLEAR) and redraw it manually. Use **Free** (Free) in **Correct/Draw** (Correct/Draw) to manually draw the margin line.

Click the **NEXT →** (NEXT) button to proceed



WORKFLOW


IMPORT

DESIGN


CAP

FILLING

Step2 Design crown



- ⑪ Check the insertion direction of the crown ,
and click the **NEXT ->** (Next) button to proceed."



- ⑫ You can change the position and shape of the crown by dragging and dropping it. You can change the type of tooth morphology from the 'Tooth Library'. You can perform adjustments from 'Adaptation to Antagonist', but please do not cut the morphology at this stage. Click the **NEXT ->** (Next) button to proceed to the next step.



WORKFLOW


IMPORT

DESIGN


CAM

MILLING



Step2 Design crown



⑬ In the anatomical tab function, you can add or remove morphology by left-dragging.



⑭ In the Free tab, you can also add or remove morphology using the  (Add/Delete) function. The  (Adapt) function allows you to adjust the occlusion with the opposing and adjacent teeth. Entering a negative value will cause interference, while entering a positive value will create a gap. After entering the values, you can perform the adjustment all at once using the  (Cut all intersections) button. The  (Smooth/Flatten) function smooths out the added or removed morphology of the tooth. Once the crown design is complete, click the  (Next) button to proceed to the next step.





WORKFLOW

IMPORT

Step3 CAM process



Step3 CAM process



15 Set up the CAM process.

The type of glass ceramic specified during patient registration will be displayed, so specify the manufacturer and size of the material to be milled."



WORKFLOW



IMPORT

DESIGN

CAM

MILLING

Step3 CAM process



⑩ If you want to change the material from the glass ceramic specified during patient registration, click the material dropdown menu from [More filters] to change the material.



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

Step3 CAM process

The sequence of screenshots illustrates the CAM process steps:

- Screenshot 1:** Shows the 'Fast' button highlighted in the 'HASS' interface.
- Screenshot 2:** Shows the 'NEXT' button highlighted in the 'Production Blank' interface.
- Screenshot 3:** Shows the 'Pin properties' dialog box with 'Tooth 45 (Anatomic crown)' selected. The 'Diameter' is set to 4. The 'Using default docking point' option is selected. The 'Avoid cavity placement' and 'Interactive mode' options are also visible.
- Screenshot 4:** Shows the 'Tools table' dialog box with a list of tools and an 'Ok' button highlighted.

3	- 0.5MM	- ZGB2-25D
4	- 1.5MM	ZGB2-75D
5	- 2.5MM	- ZGB-125D
6	- 1.0MM	- ZGB-50D

⑰ Set the machining conditions, milling position, and support size.
Click the **NEXT →** (Next) button to start the calculation.

⑱ When the tool table is displayed,
click the **Ok** (OK) button.



WORKFLOW


IMPORT

DESIGN

CAM

MILLING

Step3 CAM process



⑩ Check the “I’m done” button, and click the **NEXT →** (Next) button to finish.

*The machining files are saved in the folder: C:¥MachineTrayQueue.



WORKFLOW

IMPORT

DESIGN


CAM

MILLING

Step4 Milling



Step4 Milling



①9 Set the milling bar and material.

Import the processing data into VPanel and start the processing.



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

DWX-43W Material



Lithium Silicate

VITA SUPRINITY® PC Block PC-14
Dentsply Sirona Celtra® DUO C14 HT
Dentsply Sirona Celtra® DUO C14 LT
Lithium Silicate HT C14
Lithium Silicate LT C14

Lithium Disilicate

Cameo Glass ceramic 15.5*11*13
Cameo Glass ceramic 18*13*15
Cameo Glass ceramic 32*15*15
Cameo Glass ceramic 40*15*14
GC Initial® Lisi HT-14
GC Initial® Lisi LT-14
Lithium Disilicate HT B32/B40/B40L/C14/I12
Lithium Disilicate LT B32/C14/C16/I12
Lithium Disilicate MO C14
Lithium Disilicate MT C14
UPCERA UP-CAD Glass Ceramic HT 18*15*13
UPCERA UP-CAD Glass Ceramic HT 32*14*14
UPCERA UP-CAD Glass Ceramic HT 40*15*15

Hybrid Ceramic

SHOFU Block HC S HT/LT/S-2L
SHOFU Block HC M HT/LT/M-2L
YAMAKIN KZR-CAD HR2 S/M/L
YAMAKIN KZR-CAD HR2 GR S/L
YAMAKIN KZR-CAD HR3 M/L
VITA ENAMIC® Block EM-10
VITA ENAMIC® Block EM-14
VITA ENAMIC® multiColor EMC-14
VITA ENAMIC® multiColor EMC-16
GC CERASMART® BL 12/14/14L
GC CERASMART® HT 12/14/14L
GC CERASMART® LT 12/14/14L
VOCO Grandio® Blocs 14L HT/LT
Hybrid Ceramic HT 12/14/14L
Hybrid Ceramic LT 12/14/14L
3M Lava ULTIMATE HT 12/14L
3M Lava ULTIMATE LT 12/14L
iTENA NUMERYS HC 12/14
UPCERA Hyramic 14 HT-S/LT-S/ML



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

DWX-43W Material



Glass Ceramic

Glass Ceramic HT I10/I12/I8/V12
Glass Ceramic LT C14/I10/I12/V12
Glass Ceramic Multi C14/C14L/I12
Rosetta® BM HT C10/C12/C14
Rosetta® BM LT C10/C12/C14

Feldspar

VITA VITABLOCK® MarkII Block I12/I14/I10/I8
VITA VITABLOCK® TriLuxe forte TF-12/TF-14
VITA VITABLOCK® TriLuxe forte TF-14-14
VITA VITABLOCK® TriLuxe forte TRI-12/TRI-14
VITA VITABLOCK® TriLuxe forte TRI-14-14
VITA VITABLOCK® TriLuxe forte RL-14-14

Composite

COLTENE BRILLIANT Crios Block HT 12/14
COLTENE BRILLIANT Crios Block LT 12/14
COLTENE BRILLIANT Crios Block ST 14
ATOS BLOCK Smart Dent ATOS 14L HT
ATOS BLOCK Smart Dent ATOS 14L LT



WORKFLOW

IMPORT

DESIGN

CAM

MILLING

DWX-43W Material



AMBER MILL

Amber® Mill Direct C14
Amber® Mill HT C12/C14/C32/C40
Amber® Mill LT C12/C14/C32/C40
Amber® Mill MO C12/C14/C32/C40
Amber® Mill MT C12/C14/C32/C40



WORKFLOW

IMPORT

DESIGN

CAM

MILLING



Empowering People to Shape the Future

Make Innovation, Make Life Better

Shape Ideas / Shape the Future / Shape Change